



Comparison of Cytotoxicity between Bulk-Fill Resins and Conventional Composite Resins and the Factors Affecting Toxicity

Xinjie Lin¹, Yayi Lei¹, Jiayu Qiu^{1,2} and Jing Liu^{1*}

¹School of Stomatology, Jinan University, Guangzhou, China

²Department of Stomatology, First Affiliated Hospital of Jinan University, Guangzhou, China

Abstract

Objective: Bulk-fill resin-based composites are gradually being used in dental treatment now, but only limited data are available on their biocompatibility. The aim of this review is to analyze and compare the cytotoxicity of bulk-fill resins with that of conventional composite resins. At the same time, the mechanisms of the cytotoxic effects were analyzed from the factors of filler, monomer, and initiator.

Methods: The narrative review approach was performed. The literature search was conducted using PubMed, Web of science, and EBSCO, and non-English articles were excluded.

Results and Conclusion: The cytotoxicity of the resin was evaluated based on different biological endpoints, i.e., cell morphology, cell membrane effect, cellular metabolism, and cell growth ability. Evaluation of cytotoxicity assays can also be categorized into direct contact method, indirect contact method, and extract method according to the mode of contact. According to the results of the above cytotoxicity evaluation methods, the cytotoxicity of most bulk-fill resins at the recommended curing depth is comparable to that of conventional resins. In this review, factors such as fillers, monomers, and initiators were also discussed to analyze the mechanism of cytotoxicity.

Keywords: Bulk-fill resin; Composite resin; Cytotoxicity; Biocompatibility

OPEN ACCESS

*Correspondence:

Jing Liu, School of Stomatology, Jinan University, No.601 Huangpu Road West, Tianhe District, Guangzhou, 510630, China Tel: +86 20 85227010; E-mail: tjliu@jnu.edu.cn

Received Date: 31 Dec 2021

Accepted Date: 01 Feb 2022

Published Date: 11 Feb 2022

Citation:

Lin X, Lei Y, Qiu J, Liu J. Comparison of Cytotoxicity between Bulk-Fill Resins and Conventional Composite Resins and the Factors Affecting Toxicity. *J Dent Oral Biol.* 2022; 7(1): 1188.

ISSN: 2475-5680

Copyright © 2022 Jing Liu. This is an open access article distributed under the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

Abbreviations

HEMA: 2-Hydroxyethyl Methacrylate; TEGDMA: Triethylene Glycol Dimethacrylate; BIS-GMA: Bisphenol-A-Glycidyl-Methacrylate; BBF: Beautiful Bulk Flowable; BBR: Beautiful Bulk Restorative; hDPSCs: Human Dental Pulp Stem Cell; PRG: Pre-Reacted Glass Ionomer; NRU: Neutral Red Uptake Assay; LDH: Lactate Dehydrogenase Release Assay; OPUS: Opus Bulk Fill Flow; FBF: Filtek Bulk Fill; SDR: Smart Dentin Replacement Flow; MTT: Thiazolyl Blue Tetrazolium Bromide Assay; WST: Water-Soluble Tetrazolium Salt Assay; MTS: 3-(4,5-Dimethylthiazol-2-Yl)-5-(3-Carboxymethoxyphenyl)-2-(4-Sulfophenyl)-2H-Tetrazolium; HDPFs: Human Dental Pulp Fibroblasts; SF: Sonicfill; ROS: Reactive Oxygen Species; BisHPPP: Bis-Hydroxy-Propoxy Phenyl Propane; MA: Methacrylic Acid; UDMA: Urethane Dimethacrylate; AUDMA: Aromatic Urethane Dimethacrylate; DDDMA: 1,12-Dodecanediol Dimethacrylate; TEC: Tetric Evoceram Bulk Fill; CQ: Camphorquinone

Introduction

Light-curing composite resins are one of the most commonly used materials in dental treatment, such as cavity filling, fissure sealing, bonding, and other restorative dental treatments. However, conventional resins have substantial limitations in that they can only be placed in 2 mm increments to obtain a sufficient degree of conversion [1]. This method not only involves more operational steps but also has larger polymerization shrinkage, resulting in poorly fitting edges and microleakage [2]. Bulk-fill resin can increase the depth of cure to more than 4 mm by adding new fillers, monomers, or photoinitiators [3]. It simplifies the clinical procedure and improves the effect of restorative treatment.

In a systematic review and meta-analysis [4], follow-up results between 1 and 10 years after filling with bulk-fill resins and conventional resins showed that bulk-fill resins showed slightly lower microhardness, but they were superior to conventional resins in terms of volume shrinkage,

polymerization stress, cusp deflection, and marginal quality. Therefore, bulk-fill resins have shown similar or better performance in clinical trials compared to conventional resins and can reduce the operative time of posterior filling restorations [5], with higher patient satisfaction.

An important prerequisite for a material that can be applied safely in clinical practice is good biocompatibility. In this paper, the recent studies on the cytotoxicity of bulk-fill resins were reviewed, and the mechanisms affecting the cytotoxic effects were analyzed from the factors such as fillers, monomers, and photo initiators.

A manual search of articles in PubMed, web of science and EBSCO until January 2021 was conducted by the authors using different combinations of the terms: 'resin', 'bulk-fill', 'cytotoxicity' or 'biocompatibility'. Selection criteria included studies that evaluated the *in vitro* cytotoxicity of both conventional and bulk-filled resins simultaneously, followed by curing according to the manufacturers' instructions. A total of 369 articles were retrieved in the initial search, and only 7 articles met the inclusion criteria after full-text reading by the two authors. A comparison of bulk-fill resins and conventional resins for cytotoxicity is shown in the accompanying table.

Cytotoxicity of Bulk-Fill Resins

The main causes of cytotoxicity of composite resins are as follows:

1. Residual monomers [such as 2-Hydroxyethyl Methacrylate (HEMA), Triethylene Glycol Dimethacrylate (TEGDMA), Bisphenol-A-Glycidyl-Methacrylate (Bis-GMA), etc.]. The release of residual monomers because the degree of conversion cannot reach 100% during resins polymerization reactions.

2. Degradation products. The resin polymers will release degradation products over time. These residual monomers and degradation products can enter the oral mucosa, periodontal tissues, and dentin and affect the pulp tissue and periapical tissues through the dentinal tubule. Besides directly damaging the above-mentioned tissues, they can also cause sensitivity and pain after filling by releasing free radicals. *In vitro* cytotoxicity studies have confirmed the significant cytotoxicity of these monomers.

According to ISO 10993-5:2009 [6], the cytotoxicity evaluation of resins can be classified according to different biological endpoints or different modes of exposure.

Classification Based on Biological Endpoints

The cytotoxicity of the resin was evaluated based on different biological endpoints, i.e., cell morphology, cell membrane effect, cell metabolism ability, and cell growth ability, respectively.

Evaluation of cell morphology

The evaluation of cell morphological changes is the most intuitive method for evaluating the cytotoxicity of biomaterials, which is usually used as an auxiliary or supplement to other methods for qualitative determination. In this method, the morphological changes of cells before and after exposure to the material tested are observed by microscopy, including changes in size, shape, and nuclei, as well as the percentage of apoptotic and dead cells. Studies have shown that the reduction in cell viability and the change in cell morphology are in parallel [7,8], which means that the change in cell morphology can reflect the cytotoxicity of the material.

The cytotoxicity of most bulk-fill resins at a curing depth of 4

mm is comparable to that of conventional resins at the recommended curing depth of 2 mm [7-9], except for Beautifil bulk flowable (BBF, Shofu, Japan) and Beautifil bulk restorative (BBR, Shofu, Japan). In some studies, mouse fibroblasts [7,8] and human Dental Pulp Stem Cells (hDPSCs), were exposed to extracts of multiple cured resins for 24 h to 72 h [9]. The results showed that the cells mentioned above exposed to the extracts of BBF and BBR exhibited more cellular morphological changes of being smaller, rounder, with concentrated and fragmented nuclei, compared to other bulk-fill resins and conventional resins. This indicated that the toxicity of these two bulk-fill resins is greater. The reason is that they contain Pre-Reacted Glass ionomer (PRG) fillers [7], which release fluoride and other ions after curing. Although it can inhibit acid production by bacteria in plaque and promote enamel demineralization, fluoride has been shown to play an important role in cytotoxicity, including causing cell damage, cell cycle arrest, mitochondrial dysfunction, DNA damage, and endoplasmic reticulum stress [10]. Other bulk-fill resins without PRG fillers showed similar cytotoxicity to conventional composite resins.

Evaluation of cell membrane effects

Evaluation of cell membrane effects is another effective method to reflect the cytotoxicity of resin materials, which is reflected by the change of cell membrane permeability. Increased cell membrane permeability and compromised lysosomal membrane integrity can result in the release of some enzymes from organelles when resin materials cause cell damage. The degree of cell damage can be reflected by the Neutral Red Uptake assay (NRU assay) and the Lactate Dehydrogenase release assay (LDH assay). In particular, since the neutral red uptake assay reflects the integrity of the cellular lysosomal membrane and the toxic effect of the composite resin on lysosomes precedes the toxicity to mitochondria. The NRU assay is more sensitive than other cytotoxicity tests that reflect cellular mitochondrial damage, such as the MTT assay Nascimento et al. [11] observed by this method that the activity of cells was reduced when L929 was exposed to 11 resins for 72 h and 7 days, with the bulk-fill resin Opus Bulk Fill Flow (Opus, FGM, Brazil) and Filtek Bulk Fill (FBF, 3MESPE, Germany) showing a significant reduction in cell activity compared to the conventional resins. However, there was no significant difference in cytotoxicity between other bulk-fill resins and conventional resins. On the other hand, Haugen et al. [3] observed differences in the effects of the composite resin on different cells by LDH assay. In the study, the bulk-fill resin FBF had the greatest cytotoxicity to sensitive osteoblasts, but the toxicity of another bulk-fill resin, SureFil[®] SDR Flow (SDR, Dentsply, Germany), was slightly less than that of the conventional resin. In addition, the cytotoxic effects of these three materials on epithelial cells and fibroblasts were low and not statistically different.

Evaluation of cell metabolic activity

The degree of cell injury can be reflected by changes in the biological metabolism or the synthetic function of cells, which is the aspect of cellular metabolic activity to evaluate the cytotoxicity of resin materials. Cell metabolic activity can be detected using the method of tetrazolium salt compounds, which are degraded to colored products by the effect of mitochondrial dehydrogenases. When there are more metabolically active cells, the colored degradation products increase, so this method can be used to detect the cytotoxicity of materials. Among them, three methods, Thiazolyl Blue Tetrazolium Bromide assay (MTT assay) [11-13], Water-Soluble Tetrazolium salt assay (WST assay) [9], and 3-(4,5-dimethylthiazol-2-yl)-5-(3-

Table 1: Literature review of cytotoxicity comparison between bulk-fill resins and conventional resins.

Authors	Title	Test materials	Experimental design	Observation content	Period	Cells type	Results	Other items	Other experimental results
Toh WS, et al. [7]	<i>In Vitro</i> Biocompatibility of Contemporary Bulk-fill Composites	standard composites: ZFR/ZFF bulk-fill PRG RBCs: BBR/BBF bulk-fill non-PRG RBCs: SDR/EXP/TNC	Extracts of 2mm, 4mm thickness specimens, diluted at 1:1, 1:2, and 1:10 concentrations	① Cell Morphology ② MTS assay	24h	L929	① The cytotoxicity of PRG bulk-fill resins was greater than that of non-PRG. ② At 4-mm thickness, undiluted extracts of bulk-fill non-PRG RBCs had significantly higher cell viability than the standard composites.		
Demirel G, et al. [8]	Cytotoxic effects of contemporary bulk-fill dental composites: A real-time cell analysis	conventional flowable RBCs: CMF conventional paste-like RBCs: CME flowable bulk-fill RBCs: VBF/FBF/SDR/XTB/BBF paste-like bulk-fill RBCs: TEC/FBR/QXF/XTF/BBR	4 mm specimens in direct contact	① Cell Morphology ② Cytotoxicity assay with iCELLigence system	24, 48 and 72 h	L929	① PRG containing bulk-fill composites were severely toxic at all time points 24, 48 and 72 h. ② None of the tested composites demonstrated high cell viability (>70%) at 48 and 72 h. ③ Flowable and paste-like composites of the same brand exhibited similar cytotoxic properties.		
Lee S-M, et al. [9]	Depth-Dependent Cellular Response from Dental Bulk-Fill Resins in Human Dental Pulp Stem Cells.	conventional composite resin: ZFF bulk-fill composite resins: SDR/VBF/BBF	Extracts of 2mm, 4mm, 6mm thickness specimens	① WST assay ② Images of live and dead cells	24h	hDPSCs	All bulk-fill resins except BBF revealed cytotoxicity from 4 to 6 or 2 to 4 mm, while ZFF was cytotoxic at over 2 mm.	① Depth of Cure (scraping test and Vickers hardness profile methodology) ② Differentiation capability	① Depth of cure was detected from 3.55 to 4.02 mm in the bulk-fill resins (vs. ~2.25 mm in conventional resin), and 80% hardness compared with that of a fully polymerized top surface was determined from 4.2 to 6 mm in the bulk-fill resin (vs. 2.4 mm in conventional resin). ② Anti-differentiation was revealed at a depth of 4-6 mm in the bulk-fill resin.
Nascimento AS, et al. [11]	Physicomechanical characterization and biological evaluation of bulk-fill composite resin	conventional resins: Z350/Z350F bulk-fill resins: ABF/FBFF/FBFS/SDR/TEF/AF/XTF/XTB/Opus	extracts of the diluted resins to obtain 1:2 and 1:10 concentrations	① MTT assay ② Neutral Red test	72h and 7d	L929	① The bulk-fill resins analyzed in this study exhibited low and/or no cytotoxicity to L929 cells, except for Opus, which showed moderate cytotoxicity according to the MTT assay. ② When the NR test was used, results were not satisfactory for all composites.	① Degree of conversion (FTIR) ② Vickers hardness (SEM and AFM) ③ Surface morphology ④ Organic filler (EDS)	The resins presented acceptable values for microhardness, degree of conversion, and surface morphology.
Haugen HJ, et al. [3]	Bulk Fill Composites Have Similar Performance to Conventional Dental Composites	conventional composite material: TEC bulk fill composite materials: SDR/FBF	Extract obtained from estimated daily salivary flow rate	① Cell morphology ② LDH assay	24h	① A549 cells ② HGF ③ Primary human osteoblasts (hOB)	Conditioned media accumulated for 24 h from FBF and TEC were cytotoxic to primary human osteoblasts.	① Depth of Cure (scraping test) ② Degree of Conversion (FTIR) ③ Polymerization Shrinkage (μCT-scans) ④ Filler Mass Fraction (thermogravimetric analysis apparatus) ⑤ Morphology of Filler (SEM) ⑥ Flexural Strength and Flexural Modulus (three-point bending test) ⑦ Vickers Hardness (indentation test) ⑧ Surface Topography Parameters (abrasion test) ⑨ Monomer Elution (HPLC)	Compared to the conventional composite, the tested bulk fill materials performed equally or better in most of the tests, except for their hardness, elastic modulus, and biocompatibility with osteoblasts.
Gonçalves F, et al. [12]	A comparative study of bulk-fill composites: degree of conversion, post-gel shrinkage and cytotoxicity	conventional composites: Z350/ZFF bulk-fill composites: AB/EP/SF/FBF/FBFF/VBF	Extracts of 4mm thickness specimens	MTT assay	24h	HGF	All materials, including bulk-fill and conventional composites, were classified as non-toxic, with cell viability higher than 70%.	① Degree of conversion (FTIR) ② Volumetric shrinkage (the strain gauge method)	① Bulk-fill composites exhibited volumetric shrinkage similar to or lower (1.4 to 0.4%) than that of conventional composites (1.7%-2.1%). ② Only four of the bulk-fill composites were able to sustain a homogeneous conversion at the 4-mm depth.
Marigo L, et al. [13]	Relation between conversion degree and cytotoxicity of a flowable bulk-fill and three conventional flowable resin-composites	conventional flowable RBCs: VDF/XTE/HRI bulk fill flowable RBC: SDR	Extracts of 2 mm and 4 mm thickness specimens	MTT assay	24h	hDPSCs	All tested materials show light cytotoxic effects	① Degree of conversion (FTIR) ② Monomers Elution (HPLC) ③ Vickers Hardness (indentation test)	① DC values and the hardness change in function of thickness and type of material ② The monomers amount leached from each specimen is influenced by thickness but it is always very low.

carboxymethoxyphenyl)-2-(4-sulfophenyl)-2H-tetrazolium (MTS method) [7], showed slight cytotoxicity or no cytotoxicity for bulk-fill resins and conventional resins at the curing depth recommended by the manufacturer. However, not all bulk-fill resins could achieve proper polymerization at the recommended curing depth [11]. For example, the bulk-fill resin Opus failed to polymerize completely at 4 mm, so it had large cytotoxicity. The study also proved that the MTT assay is more sensitive than the LDH assay, but the MTT assay is unstable and the results of repeated tests are quite different.

Evaluation of cell growth ability

It mainly refers to the measurement of cell proliferation after exposure to the tested materials, and mainly includes clone formation assays. Tsuchiya et al. [14] indicated that clone formation assay is the most sensitive method for *in vitro* cytotoxicity evaluation. There is no study on the comparison of clone formation tests between bulk-filled resins and conventional resins, suggesting that this may be a direction for future research.

Classification Based on Contact Modes

Evaluations of cytotoxicity assays are categorized into direct contact method, indirect contact method, and extract method according to the mode of contact. An article [15] compared the sensitivity of these three methods and found that the extract method was poorly correlated with the other two methods and the least sensitive. However, the result of Lim et al. [16] showed the consistency of the three test models for resin cytotoxicity assay. The bulk-fill resins evaluated did not cause excessive cytotoxicity at a depth of cure of 4 mm [16].

However, the toxicity of the materials is not fully reflected by these commonly used methods, because the filling materials are separated from the pulp cells *in vivo* by dentin. The correlation between the cellular response shown by these methods and the response of the pulp cells practically *in vivo* is low. Therefore, Hume et al. [17] proposed the dentin barrier method, which can simulate the situation that the toxic substances in the material act on the pulp cells through the dentin tubules after contact with the dentin. This method has been widely used to study the cytotoxicity of dental materials in recent years. It has been shown that the cytotoxicity of zinc oxide eugenol hydromorphone is much higher than that of animal experiments when measured by existing *in vitro* cytotoxicity assays. In contrast, the cytotoxicity obtained by using the dentin barrier method was the same as that observed in clinical applications.

In addition, the cells of the traditional test methods are cultured in monolayers. Recent research [18], has established three-dimensional dentin/pulp tissue mimics as an advanced assessment tool for cytotoxicity of dental materials (DentCytoTool). It can better simulate the growth environment of cells in the oral cavity and increase the clinical relevance of the experiment.

Reasons for Differences in the Results of Cytotoxicity Evaluation

The results of the cytotoxicity tests are related to the cells used, and there are significant differences in the results of toxic effects of different cell lines on the same resin. Human Dental Pulp Fibroblasts (hDPFs) are more sensitive to the cytotoxicity tests than human Dental Pulp Stem Cells (hDPSCs) [19]. Human-derived cell lines are more sensitive to the assay of cytotoxicity than animal cell lines, such as L929 [20]. Since L929 is derived from mouse connective tissue,

which is different in species from human cells. Moreover, there are some factors such as abnormal karyotype, genetic material changes, and so on, which lead to the difference between the toxic reaction ability and that of normal cells *in vivo*. However, some studies have used human primary pulp cells [21] and primary gingival fibroblasts [22] to compare with L929 cells in toxicity tests, and the results were not significantly different. In general, a higher sensitivity indicates a more accurate response of such cells to the cytotoxicity of the material being tested. When evaluating the cytotoxicity of materials, different assays, as well as different cells, can be combined to simulate the cell growth environment *in vivo* if possible. This can reduce the experimental error and improve the correlation between the experiment and clinical treatment.

To some extent, the cytotoxicity of the resin is also related to the fluidity and viscosity of the resin. The high-viscosity type of bulk-fill resin has higher inorganic filler content and poor flow ability, but it has better mechanical properties and lower cytotoxicity after curing. In contrast, the low-viscosity type of bulk-fill resin has a lower filler content, poorer mechanical properties, and higher cytotoxicity, but has better flow ability.

According to the results of the above cytotoxicity evaluation methods, the cytotoxicity of most bulk-fill resins at the recommended curing depth is comparable to that of conventional resins, except for resins containing pre-reacted glass ionomer fillers (such as BFB, BBR, etc.) and some low-viscosity type bulk-fill resins (such as Opus, FBF, etc.), both of which are more toxic. However, the slight differences in the specific toxicity ranking may be due to the differences in the methods and evaluation criteria used in the different tests.

Factors Affecting the Cytotoxic Effect

Filler

Bulk-filled resins increase the depth of cure by reducing the filler content or increasing the filler size to reduce light scattering and increase light transmittance. Bulk-fill resin SDR [3, 23], and SonicFill (SF, Kerr, USA) [24] employ the method of increasing the filler size with irregular filler particles of about 10 μm to 20 μm in diameter to reduce light scattering and improve the degree of conversion. For the size of the filler particles, the degree of conversion is lower at wavelengths smaller than the incident blue light (400 nm to 480 nm) and closer to the output wavelength of the curing unit [25]. It is related to the scattering effect of the penetrating light in the process of activation with light. Once the filler sizes are increased over 500 nm, the degree of conversion can be improved, resulting in an increased depth of cure.

On the other hand, the filler content is negatively correlated with the light transmission [26]. For example; FBF reduces the content of filler, so as to increase the depth of cure [3]. At the same time, however, it uses different sizes of spherical high refractive index silylated zirconia/silica filler particles instead of granular glass fillers [3,24]. Although it can improve mechanical properties, such as bending strength and fracture toughness [27], the high refractive index of zirconia leads to a decrease in light transmittance [24], which makes the degree of conversion and depth of cure of FBF lower than that of SDR, and the cytotoxicity is correspondingly higher.

It is noteworthy that decreasing the filler content will bring about an increase in the toxicity of the resin. Since inorganic fillers can cause less toxicity than resin matrix [20,25], decreasing the filler content will correspondingly increase the matrix content, which leads

to a corresponding increase in unreacted monomers and degradation products. It has been proved that the reduction of filler content causes a more obvious toxic response in cells [28] and inflammation in tissues [29], which is associated with mitochondrial dysregulation caused by increased release of Reactive Oxygen Species (ROS) [30]. Materials with low filler content will release more degradation products bis-Hydroxy-Propoxy Phenyl Propane (BisHPPP, derivative of bis-GMA), Methacrylic Acid (MA), and unreacted monomer Triethylene Glycol Dimethacrylate (TEGDMA) [30] after being exposed to esterase [31]. Furthermore, TEGDMA has been proved to be closely related to inflammation, inhibition of cell proliferation and differentiation, induction of apoptosis [32], and DNA damage.

The filler content is relevant to the mechanical properties of the material [33]. Reducing the filler content will reduce its mechanical properties, such as microhardness [4] and elastic modulus [34]. Therefore, the mechanical properties, biocompatibility, and depth of cure of the materials are considered in order to obtain a suitable formulation.

Monomer

The bulk-fill resin changes the matrix composition by adding new monomers with relatively high molecular masses, such as Urethane Dimethacrylate (UDMA), Aromatic UDMA (AUDMA), and 1,12-Dodecanediol Dimethacrylate (DDDMA) [2], resulting in increased light transmission and thus increased depth of cure.

Bulk-fill resin SDR has a high degree of conversion because it contains a patented modified UDMA that chemically embeds polymerization modifiers into the resin backbone [35]. As a result, it is able to control the polymerization kinetics, improve conversion rates and reduce polymerization shrinkage [36]. Bisphenol-A-Glycidyl-Methacrylate (Bis-GMA), a matrix commonly used in conventional resins, is a highly viscous monomer that contains side hydroxyl groups and a rigid aromatic ring of bisphenol. Which negatively affects the degree of conversion [2]? While many bulk-fill resins are UDMA-based materials combining different monomers. The intermolecular hydrogen bond formed by imino groups in UDMA is weaker, which makes its viscosity lower than that of Bis-GMA. In the polymerization process, the activity and migration rate of monomers are reduced, which significantly improves the degree of conversion [37] and mechanical properties [38], without increasing polymerization shrinkage, and with less cytotoxicity [39].

According to previous studies [40], the performance of bulk-fill resins and conventional resins in leaching monomer is equivalent. However, Pongprueksa et al. [41] found that monomer release after 2 mm increment filling of conventional resins was lower than that of bulk-fill resins. This is due to the slightly lower conversion of the bulk-filled resins at 4 mm, resulting in higher monomer release than the conventional resins. Nevertheless, in a toxicity assay *in vitro*, cytotoxicity and genotoxicity of exposure to extracts were lower than exposure to single components of the composite resin [42], such as monomers, initiators, and additives, respectively. The antagonistic effect of various components in the extract [43] and the components in it can combine with saliva by protein [44], so the toxicity will be reduced. It suggests that the toxicity of the material in the actual situation is lower than that shown in the *in vitro* tests, especially the single component tests.

Other factors

Apart from changing the type of monomer and filler, the depth

of cure can be increased by adding new photoinitiators. For example, Tetric Evoceram Bulk Fill (TEC, Ivoclar-Vivadent, Liechtenstein) adds a new photoinitiator, dibenzoylgermanium derivative Ivocerin, which generates free reactive groups during the reaction, thus increasing the efficiency of the polymerization reaction [45]. Compared with traditional photoinitiators, Ivocerin has a higher degree of conversion [46,47] and color stability [48] which is considered to be an effective photoinitiator to replace Camphoroquinone (CQ) [49]. *In vitro* studies have also shown that Ivocerin has low cytotoxicity and no mutagenic effects [50,51].

In addition, the mismatch of refractive index between filler and monomer will scatter the irradiation light at the interface between resin and filler, which will lead to a lower conversion rate and greater cytotoxicity [23]. The selection of the color of the composite resin also affects the toxicity, and the composite material with a darker color has higher cytotoxicity [52].

Summary

With the rapid development of composite resin materials, more and more studies are focusing on biocompatibility while exploring the mechanical properties of the resins. In this paper, different detection methods were used to prove that the cytotoxicity of bulk-fill resins was similar to that of conventional resins, and the factors affecting the toxicity of bulk-fill resin were discussed. Although there are many reports on cytotoxicity and genotoxicity of bulk-fill resins *in vitro*, there is still a lack of research on implantation and sensitization. Therefore, future studies can perhaps concentrate on these aspects, so as to comprehensively evaluate the biocompatibility of materials and provide a guarantee for better clinical application of bulk-fill resins. Predictably, after the improvement of physical properties and mechanical properties, bulk-fill resins may be used as a substitute for conventional resins.

References

1. Abed YA, Sabry HA, Alrobeigy NA. Degree of conversion and surface hardness of bulk-fill composite versus incremental-fill composite. *Tanta Dent J*. 2015;12(2):71-80.
2. Xu T, Li X, Wang H, Zheng G, Yu G, Wang H, et al. Polymerization shrinkage kinetics and degree of conversion of resin composites. *J Oral Sci*. 2020;62(3):275-80.
3. Haugen HJ, Marovic D, Par M, Thieu MKL, Reseland JE, Johnsen GF. Bulk fill composites have similar performance to conventional dental composites. *Int J Mol Sci*. 2020;21(14):5136.
4. Cidreira Boaro LC, Pereira Lopes D, de Souza ASC, Lie Nakano E, Ayala Perez MD, Pfeifer CS, et al. Clinical performance and chemical-physical properties of bulk fill composites resin - A systematic review and meta-analysis. *Dent Mater*. 2019;35(10):e249-64.
5. Bellinaso MD, Soares FZM, Rocha RO. Do bulk-fill resins decrease the restorative time in posterior teeth? A systematic review and meta-analysis of *in vitro* studies. *J Investig Clin Dent*. 2019;10(4):e12463.
6. ISO: 10993-5. Biological evaluation of medical devices-Part 5: Tests for *in vitro* cytotoxicity. International Organization for Standardization: Geneva; 2009.
7. Toh WS, Yap AU, Lim SY. *In vitro* biocompatibility of contemporary bulk-fill composites. *Oper Dent*. 2015;40(6):644-52.
8. Demirel G, Gur G, Demirsoy FF, Altuntas EG, Yener-Ilce B, Kilicarslan MA. Cytotoxic effects of contemporary bulk-fill dental composites: A real-time cell analysis. *Dent Mater J*. 2020;39(1):101-10.
9. Lee SM, Kim SY, Kim JH, Jun SK, Kim HW, Lee JH, et al. Depth-dependent

- cellular response from dental bulk-fill resins in human dental pulp stem cells. *Stem Cells Int.* 2019;1-11.
10. Zuo H, Chen L, Kong M, Qiu L, Lu P, Wu P, et al. Toxic effects of fluoride on organisms. *Life Sci.* 2018;198:18-24.
 11. Nascimento AS, Lima DB, Fook MVL, Albuquerque MS, Lima EA, Sabino MA, et al. Physicomechanical characterization and biological evaluation of bulk-fill composite resin. *Braz Oral Res.* 2018;32:e107.
 12. Gonçalves F, Campos LMP, Rodrigues EC, Costa FV, Marques PA, Francci CE, et al. A comparative study of bulk-fill composites: Degree of conversion, post-gel shrinkage and cytotoxicity. *Braz Oral Res.* 2018;32:e17.
 13. Marigo L, Spagnuolo G, Malara F, Martorana GE, Cordaro M, Lupi A, et al. Relation between conversion degree and cytotoxicity of a flowable bulk-fill and three conventional flowable resin-composites. *Eur Rev Med Pharmacol Sci.* 2015;19(23):4469-80.
 14. Tsuchiya T, Ikarashi Y, Hata H, Toyoda K, Takahashi M, Uchima T, et al. Comparative studies of the toxicity of standard reference materials in various cytotoxicity tests and in vivo implantation tests. *J Appl Biomater.* 1993;4(2):153-6.
 15. Cao T, Saw TY, Heng BC, Liu H, Yap AU, Ng ML. Comparison of different test models for the assessment of cytotoxicity of composite resins. *J Appl Toxicol.* 2005;25(2):101-8.
 16. Lim SM, Yap A, Loo C, Ng J, Goh CY, Hong C, et al. Comparison of cytotoxicity test models for evaluating resin-based composites. *Hum Exp Toxicol.* 2017;36(4):339-48.
 17. Hume WR. A new technique for screening chemical toxicity to the pulp from dental restorative materials and procedures. *J Dent Res.* 1985;64(11):1322-5.
 18. Hadjichristou C, Papachristou E, Bonovolis I, Bakopoulou A. Three-dimensional tissue engineering-based Dentin/Pulp tissue analogue as advanced biocompatibility evaluation tool of dental restorative materials. *Dent Mater.* 2020;36(2):229-48.
 19. Tavangar MS, Attar A, Razmkhah M, Hosseini SM, Hosseini A, Monabati A, et al. Differential expression of drug resistance genes in CD146 positive dental pulp derived stem cells and CD146 negative fibroblasts. *Clin Exp Dent Res.* 2020;6(4):448-56.
 20. Schubert A, Ziegler C, Bernhard A, Burgers R, Miosge N. Cytotoxic effects to mouse and human gingival fibroblasts of a nanohybrid ormocer versus dimethacrylate-based composites. *Clin Oral Investig.* 2019;23(1):133-9.
 21. Karapinar-Kazandag M, Bayrak OF, Yalvac ME, Ersev H, Tanalp J, Sahin F, et al. Cytotoxicity of 5 endodontic sealers on L929 cell line and human dental pulp cells. *Int Endod J.* 2011;44(7):626-34.
 22. Franz A, Konig F, Skolka A, Sperr W, Bauer P, Lucas T, et al. Cytotoxicity of resin composites as a function of interface area. *Dent Mater.* 2007;23(11):1438-46.
 23. Fronza BM, Ayres A, Pacheco RR, Rueggeberg FA, Dias C, Giannini M. Characterization of inorganic filler content, mechanical properties, and light transmission of bulk-fill resin composites. *Oper Dent.* 2017;42(4):445-55.
 24. Bucuta S, Ilie N. Light transmittance and micro-mechanical properties of bulk fill vs. conventional resin based composites. *Clin Oral Investig.* 2014;18(8):1991-2000.
 25. Turssi CP, Ferracane JL, Vogel K. Filler features and their effects on wear and degree of conversion of particulate dental resin composites. *Biomaterials.* 2005;26(24):4932-7.
 26. Son SA, Park JK, Seo DG, Ko CC, Kwon YH. How light attenuation and filler content affect the microhardness and polymerization shrinkage and translucency of bulk-fill composites? *Clin Oral Investig.* 2017;21(2):559-65.
 27. Guo G, Fan Y, Zhang JF, Hagan JL, Xu X. Novel dental composites reinforced with zirconia-silica ceramic nanofibers. *Dent Mater.* 2012;28(4):360-8.
 28. Lin-Gibson S, Sung L, Forster AM, Hu H, Cheng Y, Lin NJ. Effects of filler type and content on mechanical properties of photopolymerizable composites measured across two-dimensional combinatorial arrays. *Acta Biomater.* 2009;5(6):2084-94.
 29. Olabisi Arigbiede A, Folasade Adeyemi B, Femi-Akinlosotu O. Relative biocompatibility of micro-hybrid and nano-hybrid light-activated composite resins. *J Dent Res Dent Clin Dent Prospects.* 2017;11(1):1-6.
 30. Liu B, Gan X, Zhao Y, Chen J, Yu H, Gao J, et al. TEGDMA releasing in resin composites with different filler contents and its correlation with mitochondrial mediated cytotoxicity in human gingival fibroblasts. *J Biomed Mater Res A.* 2019;107(6):1132-42.
 31. Finer Y, Santerre JP. Influence of silanated filler content on the biodegradation of bisGMA/TEGDMA dental composite resins. *J Biomed Mater Res A.* 2007;81(1):75-84.
 32. Inamitsu H, Okamoto K, Sakai E, Nishishita K, Murata H, Tsukuba T. The dental resin monomers HEMA and TEGDMA have inhibitory effects on osteoclast differentiation with low cytotoxicity. *J Appl Toxicol.* 2017;37(7):817-24.
 33. Rizzante FAP, Duque JA, Duarte MAH, Mondelli RFL, Mendonca G, Ishikiriyama SK. Polymerization shrinkage, microhardness and depth of cure of bulk fill resin composites. *Dent Mater J.* 2019;38(3):403-10.
 34. Iwasaki T, Kamiya N, Hirayama S, Tanimoto Y. Evaluation of the mechanical behavior of bulk-fill and conventional flowable resin composites using dynamic micro-indentation. *Dent Mater J.* 2021;116.
 35. Hatipoglu O, Karadas M, Er H, Turumtay EA. Effect of thermocycling on the amount of monomer released from bulk fill composite resins. *Dent Mater J.* 2019;38(6):1019-25.
 36. Ilie N, Hickel R. Investigations on a methacrylate-based flowable composite based on the SDR™ technology. *Dent Mater.* 2011;27(4):348-55.
 37. Pratap B, Gupta RK, Bhardwaj B, Nag M. Resin based restorative dental materials: Characteristics and future perspectives. *Jpn Dent Sci Rev.* 2019;55(1):126-38.
 38. Barszczewska-Rybarek IM. Structure-property relationships in dimethacrylate networks based on Bis-GMA, UDMA and TEGDMA. *Dent Mater.* 2009;25(9):1082-9.
 39. Walters NJ, Xia W, Salih V, Ashley PF, Young AM. Poly(propylene glycol) and urethane dimethacrylates improve conversion of dental composites and reveal complexity of cytocompatibility testing. *Dent Mater.* 2016;32(2):264-77.
 40. Alshali RZ, Salim NA, Sung R, Satterthwaite JD, Silikas N. Analysis of long-term monomer elution from bulk-fill and conventional resin-composites using high performance liquid chromatography. *Dent Mater.* 2015;31(12):1587-98.
 41. Pongprueksa P, De Munck J, Duca RC, Poels K, Covaci A, Hoet P, et al. Monomer elution in relation to degree of conversion for different types of composite. *J Dent.* 2015;43(12):1448-55.
 42. Yang Y, Reichl FX, Shi J, He X, Hickel R, Högg C. Cytotoxicity and DNA double-strand breaks in human gingival fibroblasts exposed to eluates of dental composites. *Dent Mater.* 2018;34(2):201-8.
 43. Wisniewska-Jarosinska M, Poplawski T, Chojnacki CJ, Pawlowska E, Krupa R, Szczepanska J, et al. Independent and combined cytotoxicity and genotoxicity of triethylene glycol dimethacrylate and urethane dimethacrylate. *Mol Biol Rep.* 2011;38(7):4603-11.
 44. Rothmund L, Shehata M, Van Landuyt KL, Schweikl H, Carell T, Geurtsen W, et al. Release and protein binding of components from resin based composites in native saliva and other extraction media. *Dent Mater.* 2015;31(5):496-504.

45. Gomes de Araujo-Neto V, Sebold M, Fernandes de Castro E, Feitosa VP, Giannini M. Evaluation of physico-mechanical properties and filler particles characterization of conventional, bulk-fill, and bioactive resin-based composites. *J Mech Behav Biomed Mater.* 2021;115:104288.
46. Alkudhairy F, AlKheraif A, Naseem M, Khan R, Vohra F. Degree of conversion and depth of cure of Ivocerin containing photo-polymerized resin luting cement in comparison to conventional luting agents. *Pak J Med Sci.* 2018;34(2):253-9.
47. Barcelos LM, Borges MG, Soares CJ, Menezes MS, Huynh V, Logan MG, et al. Effect of the photoinitiator system on the polymerization of secondary methacrylamides of systematically varied structure for dental adhesive applications. *Dent Mater.* 2020;36(3):468-77.
48. Alkudhairy F, Vohra F, Naseem M, Owais MM, Amer AHB, Almutairi KB. Color stability and degree of conversion of a novel dibenzoyl germanium derivative containing photo-polymerized resin luting cement. *J Appl Biomater Funct Mater.* 2020;18:1-7.
49. Delgado AJ, Castellanos EM, Sinhoreti M, Oliveira DC, Abdulhameed N, Geraldeli S, et al. The use of different photoinitiator systems in photopolymerizing resin cements through ceramic veneers. *Oper Dent.* 2019;44(4):396-404.
50. Moszner N, Fischer UK, Ganster B, Liska R, Rheinberger V. Benzoyl germanium derivatives as novel visible light photoinitiators for dental materials. *Dent Mater.* 2008;24(7):901-7.
51. Kowalska A, Sokolowski J, Bociog K. The photoinitiators used in resin based dental composite-A review and future perspectives. *Polymers (Basel).* 2021;13(3):470.
52. Sigusch BW, Pflaum T, Volpel A, Gretsche K, Hoy S, Watts DC, et al. Resin-composite cytotoxicity varies with shade and irradiance. *Dent Mater.* 2012;28(3):312-9.